




Snow 10 PRO

Feel More Beyond Big

Snow 10 Pro features the text-to-speech function and a 10" large screen for an excellent reading experience. Audio files and pictures can be effectively stored up to 1.5GB. A 13 mega pixels high definition camera captures high quality images and can be set self-view function providing mirror-image. A customized foldable reading stand would hold the device while reading. The device is a perfect combination of a reader, a mirror and a music player. Snow 10 Pro offers a package solution bigger than a big screen.



-  20 languages text-to-speech
-  Reader mode for easier reading
-  Self-view/distance view/close view

Key Specifications

Screen	10" touch screen; screen resolution: 1280 x 720p
HD Camera	13M pixels; auto focus/tap to focus
Color modes	Full color & 10 high contrast color modes
Viewing modes	Self/Close/Distance
Text-to-speech	Optional
Magnification	2x ~ 19x (non-text-to-speech) 1.3x ~ 8x (text-to-speech/reader mode)
Battery	Approximately 3 hours of continuous use
Fast charging time	2 hours
Port	Type-C/Micro-HDMI/Headphone
Dimensions	280 x 190 x 20.9mm (11 x 7.5 x 0.82in)
Storage	1.5GB (up to 1000 images)
Reading stand	Standard & extra foldable stand for text-to-speech
Certificates	CE FC RoHS

